

YOU ASKED 4 IT !



Clearly, it's an efficient development and optimization of fusion techniques!

Because the finest reproductibility is essential, we created the M4 Gas Fluxer for sample preparation by fusion for XRF, AA and ICP. Since 1998, it established the standard thanks to efficient agitation and superior homogenization of the melt in less time, stable flame especially at high and low temperatures, uniform heating for reproducible results, retention of volatile elements and operating without supervision. It's also easily computerized and totally user friendly. And last but not least, it's an instrument with low energy consumption and operating costs.

*Expect
the best!*

You asked for it. And you couldn't have asked for anything better!



The First and Finest in Fusion

350, rue Franquet, Sainte-Foy (Quebec) Canada G1P 4P3
Tel.: 418.656.6453 Fax: 418.656.1169 support@claisse.com
www.claisse.com

